



## FalisoMED Handicap Crank "Reha-P"

101 Reha-P, for the left leg

**Price: €224.21**

Crank Spider for rehabilitation With the rehabilitation crank attachments 101 REHA-P (for the left side) and 102 REHA P (for the right side) from FALISOMED, rehabilitation clinics or physiotherapeutic practices can respond to the limited hip or knee mobility of their patients, for example after an operation, and train their mobility again. The Reha P crank consists of a Prisma sliding rail, which is connected to the crank without movement using two screw connections and a dowel pin. The sliding block sitting on the rail holds the pedals and can be adjusted in eight stages. This allows the respective Crank Spider radius to be adjusted accordingly. The manufacturer recommends careful and precise assembly according to the instructions.

### Basic Data

Supplier Article Number:	101REHA P
EAN	999999991
MSRP:	€357.80
Fedas-Code	160362
DST-Code	2R11
Bidex-Code	202090
SB-Article	No
Brand/Supplier	FalisoMED
Unit of measure:	PIECE
Weight	0.237 kg

### Technical Data

Main colour	silver
Variant	Accessories