



## FalisoMED Handicap Crank "Reha"

102 Reha, for the right leg

**Price: 265,62 €**

Auxiliary crank with pendulum arm With the rehabilitation crank attachment Reha from FALISOMED, patients who have undergone surgery, for example, can improve their limited knee mobility by training on an Indoor Trainers or bicycle. The Crank Spider's range of motion is largely infinitely adjustable. The slide rail is available for the right (102 Reha) and left side (101 Reha) of the Crank Spider. This means that the crank arm attachment is also suitable for use in physiotherapy practices or rehabilitation clinics. The slide rail made of high-strength light metal is connected to the Crank Spider with two screw connections and a stainless steel dowel pin. A sliding block with a ball-bearing mounted pendulum crank is placed on the rail. The composition allows the pedal position to be adjusted practically infinitely and adapted to the training progress. The manufacturer emphasizes the need for careful and precise assembly according to the instructions.

### Datos maestros

Número de artículo del proveedor:	102 REHA
EAN	999999991
UVP	424,90 €
Fedas-Code	160362
DST-Code	2R11
Bidex-Code	202090
SB-Artículo	No
Marca/Proveedor	FalisoMED
Unidad de medida:	PIECE
Weight	0.355 kg

### Datos técnicos

Artículo	Accessories
Color principal	Plata